KULZER
MITSUI CHEMICALS GROUP

cara 

cara Print Clean

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 Quick Start Guide

cara makes life so easy

KURZANLEITUNG / QUICK START GUIDE

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GB Assemble Fill Place

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(GB)

LID & SHORT BUILD TABLE HANGER

- Object height: up to 170 mm
- Placing plug at lower position

LONG BUILD TABLE HANGER (INSTALL ALTERNATIVELY)

- Object height: max. 70 mm
- Placing plug at upper position

INNER CONTAINER, PLUG & MAGNET STIRRER

- Placing plug at lower position: use short table hanger, max. 2,0 liters IPA
- Placing plug at upper position: use long table hanger, max. 1,3 liters IPA
- Check bottom side foam gaskets
- Install one plug only – never both plugs in parallel
- Placing magnet stirrer at bottom of container

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PART TRAY

- If any one of the parts to be cleaned is not attached to the Build Table, the part tray must be inserted into the inner container.
- If all parts to be cleaned are attached to the Build Table, the part tray must be removed.

3 EXTRA FOAM GASKETS

- To be used as spare seals at bottom of inner container

EXTRA CONTAINER PLUG

- To be used as spare part at inner container. Install one plug only – never use both plugs in parallel

2 MAGNET STIRRER

- 1 piece to be used as spare part

cara Print Clean Quick Start Guide

Rev. 1.0 Date of release: 2019-05

Introduction: This quick start guide is intended to instruct the operator about the proper use of this device and to create awareness about the main safety aspects. This document does not substitute the Operating instructions document, which can be downloaded via the 2D barcode or the website mentioned at end of this document.

Symbols for safe use of device:

Symbol	Accompanying word(s)	Explanation
	CAUTION	This symbol indicates safety-relevant paragraphs and sections in these working instructions.
	NOTE	This symbol indicates information within the working instructions that optimizes the use of the unit.
	WARNING	This symbol warns of the dangers of electric voltage. Failure to comply with safety guidelines and procedures can result in electrical shock that may lead to injury or death.
	FLAMMABLE	This symbol warns of the dangers of flammable compounds. Failure to comply with safety guidelines and procedures can result in combustion that may lead to injury or death.
	DISCONNECT FROM MAINS	This symbol indicates danger of electric shock when unit enclosure is opened. Unplug the unit before opening it.
	READ/REFER TO INSTRUCTIONS FOR USE	This symbol directs the operator or service agent to read the Operating instructions for the Medical Equipment to find important information.
	DC POWER	This symbol indicates the Medical Equipment utilizes DC Power.
	DC JACK POLARIZED	This symbol indicates that the DC power input jack for the Medical Equipment is polarized. In particular, the DC power input jack's outer conductor is the "negative" conductor and the inner conductor is "positive" conductor.
	START BUTTON	This symbol indicates start of programmable clean cycle.
	WEEE	Disposal of old equipment according to WEEE Electrical and Electronic Equipment Act (ElektroG).

Safety information

-  **CAUTION!** Check device carefully prior to installation operation and do not work with a damaged unit under any circumstances.
-  **WARNING!** No modifications of this device are allowed.
-  **CAUTION!** Prior to operation read the cara Print Clean Operating instructions carefully. Adhere to the information provided to avoid errors, equipment damage and injury to the operator or third parties.
-  **CAUTION!** The connecting power supply cable and plug must be checked for damage prior to operation. If any damage exists, the unit must not be connected to the mains. A damaged mains power supply may only be replaced by a mains power supply of the same type. Failure to comply increases risk of electrical shock and serious harm of operator or third parties.

 **FLAMMABLE!** Isopropyl alcohol, the intended cleaning solvent, is highly flammable in both liquid and vapor states. The operating environment must be free of open flame, sparks, other ignition sources or excessive heat. Failure to comply increases risk of combustion and serious harm of operator or third parties.

 **CAUTION!** Isopropyl alcohol, the intended cleaning solvent, can be damaging to the human body through ingestion, inhalation, and bodily contact. Read and adhere to the safety information detailed on the isopropyl alcohol container's label or contact the organization responsible for your local hazardous materials regulations for more information. Mandatory safety equipment:

- Eye protection
- Suitable gloves (see chemical protective gloves, category 3, EN 374)
- Clothing that covers exposed skin
- Fire extinguisher
- Ventilation hood (air circulation in closed room/loop prohibited, free exhaust into environment only)

 **WARNING!** To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

Intended use

The cara Print Clean is an automated cleaning system that utilizes an agitated, contained volume of isopropyl alcohol (isopropanol) to clean excess uncured 3D printing liquid material from 3D printed parts ("parts"). Cleaning of parts and unit use shall be done by trained professionals in 3D printing laboratory environments. The unit is not suitable nor intended for any purpose other than cleaning in these stated industries or recreations.

Proper operation

-  **CAUTION!** The unit is suitable to be set up and operated in clean laboratory environments. The unit has been designed for continuous operation.
- Only 3D printed parts made of cured resin shall be washed in the cara Print Clean. Cleaning anything that is not a resin-based, 3D-printed part will void the warranty.
 - Do not excessively clean a 3D printed part by continuing to expose it to isopropyl alcohol past the point of removing the uncured resin. Over exposure may lead to structural degradation by weakening the part's cured resin.
 - Ensure that the outer container, inner container, inner container plug, hanging arm, Build Table, and lid are all secured in their proper positions before starting a cleaning cycle.

 **NOTE!** Isopropyl alcohol, and other officially released cleaning agents by Kulzer, are the only solvents suitable to be used in the cara Print Clean. Using any other substance will void the warranty.

 **CAUTION!** Take care to ensure that no isopropyl alcohol spills from the outer container at any time. Stray isopropylalcohol may damage the cara Print Clean's electronics or surrounding materials.

Operating instructions download

Customer is obliged to download, read and observe complete Operating instructions (www.kulzer.com):

Please find some more languages by using www.kulzer.com

© 2019 Kulzer GmbH

All Rights reserved. No parts of this work may be reproduced in any form or by any means – graphic, electronic, or mechanical, including photocopying, recording, taping, or information storage and retrieval systems – without the written permission of the publisher.

WZ1986 99001512/00 2019-05

cara makes life so easy

Find out more about our
3D Printer for dental
appliances – cara Print 4.0:
kulzer.com/cara-print

Distributor:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)