dimaprint® Print Guide

dimaprint® Print Guide

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GB Instructions for use

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SE Bruksanvisning

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KULZER
WITHE CHEMICALS GROUP

Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

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Stand: 2017-10

GB Instructions for use dima Print Guide

- 1. Introduction**
- 2. Intended use**
- 3. Contra-indication**
- 4. Hazard and Precaution**
- 5. Storage conditions, expiry date and transport**
- 6. Processing dima Print Guide**
- 7. Finishing**
- 8. Plastic and packaging waste**
- 9. Instructions for cleaning, disinfecting and sterilisation**
- 10. Delivery units**

1. Introduction
The following instructions for use are for dental technicians, dentists and oral surgeons who use dima Print Guide as a dental surgical guide material. dima Print Guide is intended exclusively for professional dental work. This instruction for use provides also information about safety and environmental aspects, a safety datasheet is available on www.kulzer.com/downloads. In case more information is needed about the processing of dima Print Guide Material contact the Kulzer Service Hotline. Also see information at the end of this document.

dima Print Guide is a monomer based on acrylic esters for manufacturing of 3D-printed surgical guides. Suitable for printing all types of surgical guides. dima Print Guide is a class I (European classification) medical device material.

dima Print Guide shall be used in combination with cara Print 4.0.

2. Intended use
dima Print Guide is a monomer mixture based on acrylic esters for manufacturing surgical guides to be used in combination with cara Print 4.0 and HiLite power 3D.

3. Contra-indication
dima Print Guide must not be used for any other purpose than dental surgical guides only. Any deviation from this instruction for use may have negative effect on the chemical and physical quality of dima Print Guide. In case of an allergic reaction, please contact a medical physician.

4. Hazard & Precaution
Inhalation: Irritating to respiratory system. High atmospheric concentrations may lead to irritation of the respiratory tract, dizziness, headache and anesthetic effects.

Skin contact: May cause sensitization by skin contact. Irritating to skin, repeated and/or prolonged contact may cause dermatitis.

Eye contact: High vapor concentration may cause irritation.

Ingestion: Low oral toxicity, but ingestion may cause irritation of the gastrointestinal tract.

Precaution: Wear protection when handling dima Print Guide. Protective glasses and nitrile gloves are advised. Information about the handling of the product can be found in the safety datasheet, which is available on www.kulzer.com/downloads.

Hazard Phrases:
H317 May cause an allergic skin reaction
H413 May cause long lasting harmful effect to aquatic life

5. Storage conditions, expiry date and transport
Store the product in the original packaging at room temperature in a dry and dark area, preferably not exceeding 25°C (77°F). Close the packaging after each use. The expiry date of the product is mentioned on the product label. In case of exceeding the expiry date, the product is no longer guaranteed in terms of treatment. Do not expose to UV-light and moisture.

6. Processing dima Print Guide
We advise to use nitrile gloves when using dima Print Guide until post-curing. Before using dima Print Guide, make sure to shake the product in the original packaging for approximately 5 minutes. Color deviations and print failures may occur when shaken insufficiently.

Make sure that the environmental conditions are as clean as possible. Dirty reservoirs or machines can cause deformation and therefore failure of the printed objects!

Pour the liquid material in the reservoir of the 3D-printing machine. Follow the instructions for use for the 3D-printer to start the printer and the printing process. When the machine has finished its program remove the building platform from machine as described in the instructions for use for the 3D-printer.

Place the platform on some paper or cloth with the built jaws facing upwards. The printed jobs can now be removed from the platform using a putty knife. Rinse the printed jobs twice with isopropanol to remove any excess material. Use an ultrasonic bath. Rinse the first time for three minutes, second time for two minutes.

The second rinse must be carried out with clean isopropanol. Rinsing with alcohol solution should not take longer than 5 minutes, as this will cause defects in the printed job. After cleaning make sure the printed part is dry and free of solvent residues. Then place the printed jobs into the HiLite power 3D for final polymerization. The final properties and also final color depend on the post-curing process. Post-curing is an UV-light treatment to ensure that dima Print materials obtain full polymer conversion, through this the residual monomer content is reduced to a minimum and the highest mechanical properties are obtained. This procedure is a necessary step to produce a biocompatible and safe end-product. We strongly advise to use the Kulzer HiLite power 3D device.

Material	Total time (min.)	Post-Curing device	Procedure
dima Print Guide	10 (2 x 5)	HiLite power 3D	Turn objects after 5 minutes

The specifications of the dima Print materials are determined on the basis of test plates, produced on our cara Print 4.0 with the corresponding settings which are registered in the software. The test plates are post-cured with the HiLite power 3D for 2 times 5 minutes wherein the objects must be turned upside down after the first 5 minute cycle. Please notice that the light sources of the light polymerisation device and the printing machine need a routine maintenance following the device instructions.

7. Finishing
Remove any support structures and finish jobs if necessary, using conventional dental methods and instruments. Differences in color nuance may occur due to

production in batches of the raw material and product or inadequate shaking of the original packaging before use or insufficient postcuring.

8. Plastic and packaging waste
The product dima Print Guide in its polymerized form is not harmful for the environment. Residual waste material in its liquid state should be delivered to a collection point for chemical waste material.

9. Instructions for cleaning, disinfecting and sterilisation

Polymerized dima 3D-printing material should be cleaned with non-chemical products. If disinfecting is required before intended use, an ethanol solution can be used. Do not use a dishwasher or thermal disinfectant. dima Print Guide can be sterilisation pouches. Make sure no mechanical forces are applied to the guide during sterilization. Let the guide cool down to room temperature before using the guide. Make sure no mechanical forces are applied to the guide during cooling down. Please make sure that the surgical guide is fully post-cured and polished before sterilization! This is important with regard to forest stability. For post-curing procedure please see "6. Processing".

10. Delivery units
The product dima Print Guide is available in the following packaging size: 1000 gr.

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