BETTER
TOGETHER

BETTER
TOGETHER

GLUMA®

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GLUMA Desensitizer & GLUMA Desensitizer PowerGel

Giving a hand to oral health.

KULZER
MITSUI CHEMICALS GROUP

Inhoud

01

Toepassing

p. 03

02

Effectiviteit

p. 06

03

Verschillen

p. 08

04

Algemeen

p. 08

05

In-vitro studies

p. 09

Toepassing

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¹ Aranha AC *et al.*: Microtensile bond strengths of composite to dentine treated with desensitizer products. J Adhes Dent. 2006, 2:85-90.

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² Mehta D, *et al.*: Efficacy of Dentin Desensitizing Agents: A Randomized Controlled Clinical Trial. *J Dent Res* 93 (Spec Iss B), 1115, 2014.

³ Schüpbach P *et al.*: Closing of dentinal tubules by Gluma desensitizer. *Eur J Oral Sci* 1997; 105: 414-421.

⁴ Qin C *et al.*: Spectroscopic investigation of the function of aqueous 2-hydroxyethylmethacrylate/glutaraldehyde solution as a dentin desensitizer. *Eur J Oral Sci* 114, 2006:354-9.

⁵ Ishihata H *et al.*: In vitro dentin permeability after application of Gluma® desensitizer as aqueous solution or aqueous fumed silica dispersion. *J Appl Oral Sci* 19(2), 2011:147-53.

⁶ Sabatini C *et al.*: Inhibition of endogenous human dentin MMPs by Gluma. *Dental Mat* 30, 2014: 752-8.

⁷ Felton D *et al.*: Inhibition of bacterial growth under composite restorations following GLUMA pretreatment. *JDR*, 68 (3), 1989: 491-5.

⁸ Bedran-Russo AK *et al.*: Changes in stiffness of demineralized dentin following application of collagen cross-linkers. *J of Biomedical Materials Research Part B: Applied Biomaterials*, 86 (B), 2008: 330-4.

⁹ Lopes AO *et al.*: Evaluation of different treatment protocols for dentine hypersensitivity: an 18-month randomized clinical trial. *Lasers Med Sci*, 32, 2017:1023-30.

¹⁰ Lopes AO *et al.*: Evaluation of different treatment protocols for dentine hypersensitivity: an 18-month randomized clinical trial. *Lasers Med Sci*, 32, 2017:1023-30.

¹¹ Silva SMA *et al.*: Effect of Oxalate Desensitizer on the Durability of Resin-Bonded Interfaces. *Operative Dentistry* 35-6, 2010: 610-617.

¹² Acar O *et al.*: The effect of dentin desensitizers and Nd:YAG laser pre-treatment on microtensile bond strength of self-adhesive resin cement to dentin. *J Adv Prosthodont* 6, 2014: 88-95.

¹³ Dondi Dall' Orologio *et al.*: In vitro and in vivo evaluation of the effectiveness of three dentin desensitizing treatment regimens. *American Journal of Dentistry* 27 (3), 2014: 139-144.

¹⁴ Mehta D *et al.*: Randomized controlled clinical trial on the efficacy of dentin desensitizing agents. *Acta Odontologica Scandinavica*. 2014; Early Online, 1-6.

¹⁵ Vora J *et al.*: Effects of two topic desensitizing agents and placebo on dentin hypersensitivity. *AJD* 25,5, 2012:293-8.

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In-vitro studies

- 1 **Sivaramakrishnan G, Sridharan K:** Fluoride varnish versus glutaraldehyde for hypersensitive teeth: a randomized controlled trial, meta-analysis and trial sequential analysis. *Clin Oral Investig.* 2018 Apr 2. doi: 10.1007/s00784-018-2428-8. [Epub ahead of print]
- 2 **Diniz A, Lima S, Tavarez R, Borges AH, Pinto S, Tonetto MR, Loguercio AD, Bandéca MC:** Preventive Use of a Resin-based Desensitizer Containing Glutaraldehyde on Tooth Sensitivity Caused by In-office Bleaching: A Randomized, Single-blind Clinical Trial. *Oper Dent.* 2018 Mar 23. doi: 10.2341/17-020-C. [Epub ahead of print]
- 3 **Hajizadeh H, Nemati-Karimooy A, Majidinia S, Moeintaghavi A, Ghavamnasiri M:** Comparing the effect of a desensitizing material and a self-etch adhesive on dentin sensitivity after periodontal surgery: a randomized clinical trial. *Restor Dent Endod.* 2017 Aug;42(3):168-175. doi: 10.5395/rde.2017.42.3.168. Epub 2017 Jul 21.
- 4 **Idon PI, Esan TA, Bamise CT:** Efficacy of Three In-Office Dentin Hypersensitivity Treatments. *Oral Health Prev Dent.* 2017;15(3):207-214. doi: 10.3290/j.ohpd.a38523.
- 5 **Lopes AO, de Paula Eduardo C, Aranha ACC:** Evaluation of different treatment protocols for dentin hypersensitivity: an 18-month randomized clinical trial. *Lasers Med Sci.* 2017 Jul;32(5):1023-1030.
- 6 **Kara HB, Cakan U, Yilmaz B, Inan Kurugol P:** Efficacy of Diode Laser and Gluma on Post-Preparation Sensitivity: A Randomized Split-Mouth Clinical Study. *J Esthet Restor Dent.* 2016 Nov 12;28(6):405-411.
- 7 **Samuel SR, Khatri SG, Acharya S, Patil ST:** Evaluation of instant desensitization after a single topical application over 30 days: a randomized trial. *Aust Dent J.* 2015 Sep;60(3):336-42. doi: 10.1111/adj.12341. Epub 2015 Jul 24.
- 8 **Patil SA, Naik BD, Suma R:** Evaluation of three different agents for in-office treatment of dentinal hypersensitivity: a controlled clinical study. *Indian J Dent Res.* 2015 Jan-Feb;26(1):38-42.
- 9 **Samuel SR, Khatri SG, Acharya S:** Clinical Evaluation of self and professionally applied desensitizing agents in relieving dentin hypersensitivity after a single topical application: A Randomized Controlled Trial. *J Clin Exp Dent.* 2014 Oct 1;6(4):e339-43.
- 10 **Dall'Orologio GD, Ishihata H, Finger WJ, Sasaki K:** In vitro and in vivo evaluation of the effectiveness of three dentin desensitizing treatment regimens. *Am J Dent.* 2014 Jun;27(3):139-44.
- 11 **Ding YJ, Yao H, Wang GH, Song H:** A randomized double-blind placebo-controlled study of the efficacy of Clinpro XT varnish and Gluma dentin desensitizer on dentin hypersensitivity. *Am J Dent.* 2014 Apr;27(2):79-83.

- 12 Mehta D, Gowda VS, Santosh A, Finger WJ, Sasaki K: Randomized controlled clinical trial on the efficacy of dentin desensitizing agents. *Acta Odontol Scand.* 2014 Nov;72(8):936-41.
- 13 Femiano F, Femiano R, Lanza A, Festa MV, Rullo R, Perillo L: Efficacy of diode laser in association to sodium fluoride vs Gluma desensitizer on treatment of cervical dentin hypersensitivity. A double blind controlled trial. *Am J Dent.* 2013 Aug;26(4):214-8.
- 14 Lopes AO, Eduardo Cde P, Aranha AC: Clinical evaluation of low-power laser and a desensitizing agent on dentin hypersensitivity. *Lasers Med Sci.* 2015 Feb;30(2):823-9.
- 15 Mehta D, Venkata S, Naganath M, LingaReddy U, Ishihata H, Finger WJ: Clinical trial of tooth desensitization prior to in-office bleaching. *Eur J Oral Sci.* 2013 Oct;121(5):477-81.
- 16 Lopes AO, Aranha AC: Comparative evaluation of the effects of Nd:YAG laser and a desensitizer agent on the treatment of dentin hypersensitivity: a clinical study. *Photomed Laser Surg.* 2013 Mar;31(3):132-8.
- 17 Vora J, Mehta D, Meena N, Sushma G, Finger WJ, Kanehira M: Effects of two topical desensitizing agents and placebo on dentin hypersensitivity. *Am J Dent.* 2012 Oct;25(5):293-8.
- 18 Ehlers V, Ernst CP, Reich M, Kämmerer P, Willershausen B: Clinical comparison of gluma and Er:YAG laser treatment of cervically exposed hypersensitive dentin. *Am J Dent.* 2012 Jun;25(3):131-5.
- 19 Guentsch A, Seidler K, Nietzsche S, Hefti AF, Preshaw PM, Watts DC, Jandt KD, Sigusch BW: Biomimetic mineralization: long-term observations in patients with dentin sensitivity. *Dent Mater.* 2012 Apr;28(4):457-64.
- 20 Brahmabhatt N, Bhavsar N, Sahayata V, Acharya A, Kshatriya P: A double blind controlled trial comparing three treatment modalities for dentin hypersensitivity. *Med Oral Patol Oral Cir Bucal.* 2012 May 1;17(3):e483-90.
- 21 Sethna GD, Prabhuji ML, Karthikeyan BV: Comparison of two different forms of varnishes in the treatment of dentine hypersensitivity: a subject-blind randomised clinical study. *Oral Health Prev Dent.* 2011;9(2):143-50.
- 22 Yu X, Liang B, Jin X, Fu B, Hannig M: Comparative in vivo study on the desensitizing efficacy of dentin desensitizers and one-bottle self-etching adhesives. *Oper Dent.* 2010 May-Jun;35(3):279-86.
- 23 Aranha AC, Pimenta LA, Marchi GM: Clinical evaluation of desensitizing treatments for cervical dentin hypersensitivity. *Braz Oral Res.* 2009 Jul-Sep;23(3):333-9.

- 24 **Ozen T, Orhan K, Avsever H, Tunca YM, Ulker AE, Akyol M:** Dentin hypersensitivity: a randomized clinical comparison of three different agents in a short-term treatment period. *Oper Dent.* 2009 Jul-Aug;34(4):392-8.
- 25 **Jalalian E, Meraji N, Mirzaei M:** A comparison of the efficacy of potassium nitrate and Gluma desensitizer in the reduction of hypersensitivity in teeth with full-crown preparations. *J Contemp Dent Pract.* 2009 Jan 1;10(1):66-73.
- 26 **de Assis Cde A, Antoniazzi RP, Zanatta FB, Rösing CK:** Efficacy of Gluma Desensitizer on dentin hypersensitivity in periodontally treated patients. *Braz Oral Res.* 2006 Jul-Sep;20(3):252-6.
- 27 **Kakaboura A, Rahiotis C, Thomaidis S, Doukoudakis S:** Clinical effectiveness of two agents on the treatment of tooth cervical hypersensitivity. *Am J Dent.* 2005 Aug;18(4):291-5.
- 28 **Sobral MA, Garone-Netto N, Luz MA, Santos AP:** Prevention of postoperative tooth sensitivity: a preliminary clinical trial. *J Oral Rehabil.* 2005 Sep;32(9):661-8.
- 29 **Duran I, Sengun A:** The long-term effectiveness of five current desensitizing products on cervical dentine sensitivity. *J Oral Rehabil.* 2004 Apr;31(4):351-6.
- 30 **Dondi dall'Orologio G, Lorenzi R, Anselmi M, Opisso V:** Dentin desensitizing effects of Gluma Alternate, Health-Dent Desensitizer and Scotchbond Multi-Purpose. *Am J Dent.* 1999 Jun;12(3):103-6.
- 31 **Dondi dall'Orologio G, Malferrari S:** Desensitizing effects of Gluma and Gluma 2000 on hypersensitive dentin. *Am J Dent.* 1993 Dec;6(6):283-6.
- 32 **Schüpbach P, Lutz F, Finger WJ:** Closing of dentinal tubules by Gluma desensitizer. *Eur J Oral Sci.* 1997 Oct;105(5 Pt 1):414-21.
- 33 **Quarnstrom F, Collier N, McDade E, McLean K, Munk A, Nicholls J:** A randomized clinical trial of agents to reduce sensitivity after crown cementation. *Gen Dent.* 1998 Jan-Feb;46(1):68-74
- 34 **Sivaramkrishnan G, Sridharan K:** Fluoride varnish versus glutaraldehyde for hypersensitive teeth: a randomized controlled trial, meta-analysis and trial sequential analysis. *Clin Oral Investig.* 2018 Apr 2. doi: 10.1007/s00784-018-2428-8. [Epub ahead of print]

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Kulzer Benelux B.V.