
SOPIRA® Carpule® Syringe Non-Aspirating

- DE** Gebrauchsanweisung
- GB** Instructions for use
- FR** Mode d'emploi
- ES** Instrucciones de uso
- IT** Istruzioni per l'uso
- NL** Gebruiksaanwijzing
- SE** Bruksanvisning
- DK** Brugervejledning
- NO** Bruksanvisning
- FI** Käyttöohjeet
- EE** Kasutusjuhend

Manufacturer:
Henke-Sass, Wolf GmbH
Kaltenstraße 1
78532 Tuttlingen (Germany)

Distributor:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)
Made in Germany

MD
Medizinprodukt
Medical device
Dispositif médical
Producto sanitario
Dispositivo médico
Medisch hulpmiddel

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GB SOPIRA Carpule Syringe
Non-aspirating dental ampoule syringe for 1.8 ml ampoules with metric threads. The syringe is intended for applying anaesthetic in cylindrical ampoules.
Range of indications: Periodontology, endodontics and restorative dentistry.
All standard needles can be used with the screw-on front section (1).

Contra-indications
The use of this product is contraindicated in case of known or suspected allergies against components of this product.
Side effects
This product or one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.
Special instructions: Only to be used by dentists and for its intended use.

Handling the syringe
1) Fully retract the plunger rod (5) against the force of the plunger spring and rotate the syringe barrel (2) 90°. (Fig. B1)
2) Insert the cylindrical ampoule. (Fig. B2)
3) Retract the plunger rod (5) against the force of the plunger spring and return the syringe barrel to its normal position.
4) Screw the needle into the screw-on front section (1). The end of the needle pierces the rubber membrane of the cylindrical ampoule. The syringe is now ready to use.

Preparation
The syringe must be prepared prior to initial use and each subsequent use (syringes are supplied clean but not sterile).

Preparation at the place of use:	The syringe must be checked for any damage such as coating each time before use. Syringes must not longer be used if they are damaged. Remove debris from the instruments immediately after use. Do not use any fixing agents or hot water (>40°C) as this will fix residues and affect the cleaning results.
Transport:	The medical devices should be stored securely and transported in a closed container to avoid damage and prevent environmental contamination.
Manual pre-cleaning:	Clean the medical devices with a soft multi-purpose brush under cold tap water until all visible residues and debris have been removed. Rise lumens, holes and threads with a water jet for at least 10 seconds at a pressure of 3.8 bar

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Cleaning:	Place the medical devices on a sieve/tray on the trolley/cart and start the cleaning process. 1. Prewash for 4 mins with cold water 2. Drain 3. Pre-wash for 10 mins with 0.5% neodisher MediClean, Dr. Weigert (Hamburg) 4. Drain 5. Neutralize for 6 mins with 0.1% neodisher Z, Dr. Weigert (Hamburg) (> 40°C) 6. Drain 7. Rinse for 3 mins with deionized water (<40°C) 8. Drain
Disinfection:	Perform automatic thermal disinfection in compliance with national AD value requirements (see ISO 15883).
Drying:	Use the drying cycle on the cleaning/disinfection device to dry the exterior of the medical devices A hot-air cloth can be used for additional manual drying. Cavities in medical devices must be dried with sterile compressed air.
Functional testing and maintenance:	Perform a visual inspection for cleanliness, followed by maintenance and functional testing in accordance with the operating instructions. When necessary, repeat the preparation process until the medical device is visually clean.
Packaging:	Medical devices must be packaged for sterilization in accordance with ISO 11607 and EN 868
Sterilization (autoclaving):	Sterilize the devices using the fractionated pre-vacuum method (in accordance with ISO 13060/ISO 17665) in compliance with national requirements. 1. Fractionated pre-vacuum (3x) 2. Sterilization temperature of 134°C 3. Minimum holding time: 3 mins (full cycle) 4. Drying time: at least 10 mins
Storage:	Sterilized medical devices must be stored in a dry, clean and dust-free environment at moderate temperatures between +5°C and +40°C.
Validating the preparation process:	The following materials and machines were used in the validation process: Detergent: neodisher MediClean (alkaline); Dr. Weigert; Hamburg Neutralising agent: neodisher Z; Dr. Weigert; Hamburg Cleaning/Disinfecting device: Miele G 7836 CD

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Stand: 2020-04

Additional instructions:	The German Medical Devices Operator Ordinance (MPBetreibV) stipulates that users are responsible for validating the preparation processes they undertake. This applies even if the chemicals and machines described above are not available. The user must ensure that the preparation process, including resources, materials and personnel, is suitable for achieving the required results. Industry standards and national legislation require the use of validated processes.
---------------------------------	---

Do not loosen or remove any screws (pivot mechanism/finger ring) before starting the preparation process.

According to the EU Medical Devices Regulation, users / patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred. For disposal please follow the instructions on the material safety data sheet and national regulations. Please quote batch number and the article number in all correspondence about the product.

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