



cara I-Butment® regular/angled
Platform overview/Plattformübersicht

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

cara I-Butment® regular/angled

Implant System compatible Implantatsystem kompatibel	Specification Kennzeich- nung	Platform Plattform	cara I-Butment regular			cara I-Butment angled			
			Ti	Zr/int. (two-piece/ zweiteilig)	CoCr	Ti	CoCr	max. angulation	Torque Ncm
Astra Tech EV	EV-Series	Ø3.0 mm	●	-	-	-	-	-	25
		Ø3.6 mm	●	-	-	-	-	-	
		Ø4.2 mm	●	-	-	-	-	-	
		Ø4.8 mm	●	-	-	-	-	-	
		Ø5.4 mm	●	-	-	-	-	-	
Astra Tech OsseoSpeed®	S-Series	Ø3.5/4.0 mm	●	●	●	●	●	15°	20
		Ø4.5/5.0 mm	●	●	●	●	●	5°	
Bego® Semados®	B-Series	Ø3,25 mm	●	-	-	-	-	-	30
		Ø3,75 mm	●	-	-	-	-	-	
		Ø4,10 mm	●	-	-	-	-	-	
		Ø4,50 mm	●	-	-	-	-	-	
		Ø5,50 mm	●	-	-	-	-	-	
Biomet 3i™ Certain®	H-Series	Ø3.4 mm	●	●	●	●	●	15°	20
		Ø4.1 mm	●	●	●	●	●	15°	
		Ø5.0 mm	●	●	●	●	●	15°	
Biomet 3i™ Osseotite®	I-Series	Ø3.4 mm	●	●	●	●	●	20°	35
		Ø4.1 mm	●	●	●	●	●	20°	
		Ø5.0 mm	●	●	●	●	●	20°	
Brånemark System®	K-Series	NP	●	●	●	●	●	20°	35
		RP	●	●	●	●	●	25°	
		WP	●	●	●	N/A	N/A	N/A	
Bredent Sky®	Sky- Series	Ø3,5 mm	●	-	-	-	-	-	25
		Ø4,0 mm	●	-	-	-	-	-	
		Ø4,5 mm	●	-	-	-	-	-	
Camlog®	C-Series	Ø3.3 mm	●	●	-	-	-	-	20
		Ø3.8 mm	●	●	-	-	-	-	
		Ø4.3 mm	●	●	-	-	-	-	
		Ø5.0 mm	●	●	-	-	-	-	
		Ø6.0 mm	●	●	-	-	-	-	
Conelog®	-	Ø3.3 mm	●	-	-	-	-	-	20
		Ø3.8/4.3 mm	●	-	-	-	-	-	
		Ø5.0 mm	●	-	-	-	-	-	
Medentis® ICX	-	Ø3.45 – 4.8 mm	●	-	-	-	-	30	
NEOSS®	-	Ø3.5 – 5.5 mm	●	-	-	-	-	32	
NobelActive®	F-Series	NP Ø3.5 mm	●	●	●	●	●	15°	35
		RP Ø4.3/5.0 mm	●	●	●	●	●	10°	

cara I-Butment® regular/angled

Implant System compatible Implantatsystem kompatibel	Specification Kennzeichnung	Platform Plattform	cara I-Butment regular			cara I-Butment angled			
			Ti	Zr/int. (two-piece/ zweiteilig)	CoCr	Ti	CoCr	max. angulation	Torque Ncm
NobelReplace®	E-Series	NP	•	•	•	•	•	5°	35
		RP	•	•	•	•	•	5°	
		WP	•	•	•	•	•	5°	
		Ø6.0 mm	•	•	•	•	•	5°	
Straumann® Bone Level	L-Series	NC Ø3.3 mm	•	•	•	•	•	15°	35
		RC Ø4.1/4.8 mm	•	•	•	•	•	15°	
Straumann® SynOcta®	N-Series	NN Ø3.3 mm	•	•	•	•	•	15°	35
		RN Ø4.8 mm	•	•	•	•	•	15°	
		WN Ø6.5 mm	•	•	•	•	•	15°	
Sweden & Martina®	SM-Series	Ø3.3 mm	•	-	-	-	-	-	25
		Ø3.8 mm	•	-	-	-	-	-	
		Ø4.25 mm	•	-	-	-	-	-	
		Ø5.0 mm	•	-	-	-	-	-	
Thommen® Implant System	-	Ø3.5 mm	•	•	•	•	•	15°	15
		Ø4.0 mm	•	•	•	•	•	15°	20
		Ø4.5 mm	•	•	•	•	•	15°	
		Ø5.0 mm	•	•	•	•	•	15°	
		Ø6.0 mm	•	•	•	•	•	15°	25
Dentsply Implants Xive®/Frialit	T-Series	Ø3.4 mm	•	•	•	•	•	10°	25
		Ø3.8 mm	•	•	•	•	•	10°	
		Ø4.5 mm	•	•	•	•	•	10°	
		Ø5.5 mm	•	•	•	•	•	10°	
Zimmer Dental TaperedScrew-Vent®	R-Series	Ø3.5 mm	•	•	•	•	•	15°	30
		Ø4.5 mm	•	•	•	•	•	15°	
		Ø5.7 mm	•	•	•	•	•	15°	



cara I-Butment is available in Zr to the following platforms:

- Regular/Angled: Brånemark system Implant level NP/RP.
- Regular: Brånemark system Implant level WP.

All products are not available in all markets. Nicht alle Produkte sind in allen Ländern verfügbar.

Brånemark System®, NobelReplace® and NobelActive® are trademarks of the Nobel Biocare group. Straumann® and Synocta® are trademark of the Straumann Group. Astra Tech OsseoSpeed®, Xive®, Frialit are trademarks of the Dentsply Group. 3i®, Biomet®, Biomet 3i™, Certain®, Osseotite® are trademarks of Biomet 3i, Inc. Screw-Vent® is a trademark of Zimmer Dental, USA. Thommen® Implant System are trademarks of Thommen Medical Group. Camlog® and Conelog® are trademarks of Camlog Biotechnologies AG. Medentis® ICX-templant® is a trademark of medendis medical GmbH. NEOSS® is a trademark of NEOSS Ltd. Bego® Semados® is a trademark of Bego Bremer Goldschlägerei GmbH & Co. KG. Sky® is a trademark of Bredent medical GmbH & Co. KG. Sweden&Martina® is a trademark of Sweden&Martina spa.

All trademarks and company names are the property of their respective owners. Die eingetragenen Warenzeichen und Firmennamen sind Eigentum der jeweiligen Firma.

Our specialists in the field will be pleased to provide you with information and service. Please contact our local Sales Representative in your country.

ORDER & CUSTOMER SERVICE

☎ CARA-SERVICE@KULZER-DENTAL.COM

BESTELL- UND SERVICE-HOTLINE

(DEUTSCHLAND):

☎ **0800.437 25 22**

(gebührenfrei)



Contact in Germany

Kulzer GmbH
Leipziger Straße 2
63450 Hanau, Germany
cara-service@kulzer-dental.com