










dima Print Materials

Light-curing resins tailored to the needs of each indication.

dima Print materials ensure excellent clinical results and high efficiency for each indication. The cara print system and the materials are perfectly geared to each other.

	Model		Gingiva mask	Nightguard/ splint	Casting structure	Individual impression tray	Implant drilling guide
							
Product name	dima Print Stone beige	dima Print Stone teal	dima Print Gingiva Mask	dima Print Splint clear	dima Print Cast ruby	dima Print Impression blue	dima Print Guide
Indication	<ul style="list-style-type: none"> Situation model Section model Ortho model Implant model 	<ul style="list-style-type: none"> Orthodontic model 	<ul style="list-style-type: none"> Gingiva mask for models 	<ul style="list-style-type: none"> Functional splint Nightguard 	<ul style="list-style-type: none"> C&B (up to 5 units) Partial denture Press ceramic 	<ul style="list-style-type: none"> Impression tray 	<ul style="list-style-type: none"> Surgical dental drilling guides
Features	<ul style="list-style-type: none"> High precision for implant models Heat resistance suitable for deep drawing Aesthetic color for veneering Matte surface finish Color: gypsum beige 	<ul style="list-style-type: none"> Extra high heat resistance, ideal for deep drawing Strong contrast to metal wires High surface hardness Extra high printing speed Easy to grind if necessary Color: teal 	<ul style="list-style-type: none"> Reproduction of fine details Elastic deformation/behavior Flexible Tear-proof Silicon-free Color: pink 	<ul style="list-style-type: none"> High accuracy High elongation at break High wearing comfort Color: transparent clear 	<ul style="list-style-type: none"> Very fluid with no sedimentation Flexible clasps Highly precise and clean crown edges Residue-free burnout for good casting results Comparable to milled cast material (PMMA) Smooth metal/cast surfaces Color: red 	<ul style="list-style-type: none"> High layer thickness and fast printing Color: blue 	<ul style="list-style-type: none"> High precision Indicator for sterilization by changing color High stiffness Color: orange, after sterilization light transparent orange
MD class	Not applicable	Not applicable	Not applicable	Ila	Not applicable	I	I
							
Product name				dima Print Ortho	dima Print Cast		
Indication				<ul style="list-style-type: none"> Splint 	<ul style="list-style-type: none"> Partial denture 		
Features				<ul style="list-style-type: none"> High accuracy High flexural strength Color: transparent blue 	<ul style="list-style-type: none"> High precision High green stability Burn out without residue Color: raspberry 		
MD class				Ila	Not applicable		